



Flake search station user's manual

Note: Do not rename or delete the zlib123dllx64 folder in desktop/Tom since it's part of the onnxruntime-gpu environment

1. Software Set up:

- X20 Magnification
- Turn the NIS-Element imaging software on
- Adjust brightness level and lighting
- Focus the image
- Press “Auto White”
- Make sure the image is fully zoomed in with no black edges




- Upload the corresponding script to the scanning stage, and disconnect it.
- I recommend to set the exposure time to less than 3ms to avoid motion blur

Flake search station user's manual

How to upload:

- Select the dimension you want
- Double click to open the INO file
- Specify the USB Port here
- You may have to specify the model and USB port manually. The board is Arduino UNO, and the port is the port with a (USB) after the name
- Upload it by pressing this button
- A “Done Uploading” notification will pop up after a few seconds
- Disconnect the board from power after uploading




2. Scanning Stage Setup

- Connect the scanning stage to the microscope with rubber band.

Flake search station user's manual

- Make sure the rubber bands have enough free space in the direction of the scan (knots should not touch the motors or the microscope stage. Note for future upgrade: Replace the connection material or “weld” them together so no knots are used)





- Put 1-2 plastic “spacers” to ensure tension, the end of the platform should be slightly in the air.
- Focus the substrate, note the start point of the scan (I recommend to start with the top right corner)



Flake search station user's manual

3. Scanning Script

- Double click the “Flake Detector” script in the folder.
- Make sure the NIS-Elements software is at the foreground with nothing blocking it
- A “Confirm if ready” window will show, Press “Ok” once all previous settings are complete.
- After the user interface shows up (takes ~20 seconds to start). Make sure it doesn't block the software. Specify the region of interest (default value is 8mmx8mm) by clicking the button. A message will pop up



to confirm you've set the dimension

- Press “Start Scanning” on the keyboard while connecting the scanning stage to power.

Flake search station user's manual

- Scan will start automatically the second the stage is connect to the PC
- Press “Stop Scanning” after the scan stops. Report and images are saved automatically
- Close the py.exe window after it completes the scan.

4. Output

- Script will stop running after pressing “Stop scanning” button
- Saved Flake Images are in the “Flakes” folder
- Manually Saved Images are in the directory folder
- Image Coordinates are recorded in the “report.txt” text file in the directory folder.
- Coordinates are with respect to the start point, the x axis points leftwards, the y axis points downwards.